

物理实验教学中心

Physics Experiment Center



Torsional pendulum


Li Bin
NJUPT

Experimental Purposes:

1. Determination of the **torsional constant** K of spiral spring using torsion pendulum;
2. Determination of the **rotational inertia** of different objects with a torsion pendulum, and compared with the theoretical value;
3. Verification of the **parallel axis theorem**.



Structure of torsional pendulum



Principles

1. Rotational inertia: I

2. Torsional pendulum

When the object is turned in a horizontal plane at a certain angle θ , the spring restoring moment M :

$$M = -K \theta \quad K \text{ is the torsional constant .}$$


And we also know that : $M = I \cdot \beta$,

where $\beta = d^2 \theta / dt^2 = -K \theta / I = -\omega^2 \theta$,


the vibration period $T = 2\pi / \omega = 2\pi (I/K)^{1/2}$

3. Parallel axis theorem

$$I_c = I_0 + mx^2$$



Instruments:



A hollow metal cylinder, a solid plastic cylinder, a metal rod, and two small metal cylinders

Instrument panel



Adjustment and use:

1. Horizontal adjustment using air-bubble level;

2. Baffle rod can cover the infrared emitting and receiving hole of the photoelectric detector.


3. Swing angle should be between 40 - 90 degrees.

Steps:


1. Determine the vibration periods ($10T$) of the loading plate and with cylindrical, calculate the torsion constant K ;



2. Determine the vibration periods ($10T$) of the hollow metal cylinder and the mental rod.



3. Verification of the parallel axis theorem;



4. Other measurements:

Electronic balance-----mass,

Tape-----length of the metal rod,

Vernier caliper-----Inside and outside diameters.

Post-processing:

1. Calculate K;

2. Calculate the experimental and theoretical rotational inertia I_{exp} and I_{theo} ,

3. Calculate the percentage error η :

$$\eta = |I_{\text{exp}} - I_{\text{theo}}| / I_{\text{theo}} * 100\%$$

Rotational inertia: $k = 4\pi^2 \frac{I_1'}{\bar{T}_1^2 - \bar{T}_0^2} = \text{_____ kg} \cdot \text{m}^2 \cdot \text{s}^{-2}$

Mass of small metal cylinder: $m = \text{_____ g}$

Geometry of small metal cylinder:

$D_o = \text{_____ mm}$, $D_i = \text{_____ mm}$, Height $L_1 = \text{_____ mm}$.

Table I. Rotational inertia

Object	Mass/kg	Geometry	Period	$I_{\text{theo}} / (\text{kg} \cdot \text{m}^2)$	$I_{\text{exp}} / (\text{kg} \cdot \text{m}^2)$	Percentage error
Plate	/	/	$10T/s$	/	$\left(I_0 = \frac{I_1' \bar{T}_0^2}{\bar{T}_1^2 - \bar{T}_0^2} \right)$	/
			\bar{T}_0/s			
Plastic cylinder		D_1/mm	$10T/s$	$\left(I_1' = \frac{1}{8} m D_1^2 \right)$	$\left(I_1 = \frac{K \bar{T}_1^2}{4\pi^2} - I_0 \right)$	
		D_1/mm	\bar{T}_1/s			
Hollow cylinder		D_{out}/mm	$10T_2/s$	$\left[I_2' = \frac{1}{8} m (\bar{D}_{\text{out}}^2 + \bar{D}_{\text{in}}^2) \right]$	$\left(I_2 = \frac{K \bar{T}_2^2}{4\pi^2} - I_0 \right)$	
		D_{out}/mm				
		D_{in}/mm				
		D_{in}/mm	\bar{T}_2/s			
Metal rod		L/mm	$10T_4/s$	$\left(I_4' = \frac{1}{12} m L^2 \right)$	$\left(I_4 = \frac{K}{4\pi^2} \bar{T}_4^2 \right)$	
			\bar{T}_4/s			

Table II. Verification of the parallel axis theorem

Position $x/10^{-2}$ m	5.00	10.00	15.00	20.00	25.00
10 periods $10T/s$					
Period \bar{T}/s					
Exp. value/kg \cdot m ² $I = \frac{k}{4\pi^2} \bar{T}^2$					
Theo. value/kg \cdot m ² $I' = I'_4 + 2mx^2 + I'_5$					
percentage error					

$$I'_5 = 2 \left[\frac{1}{16} m(D_0^2 + D_1^2) + \frac{1}{12} mL_1^2 \right]$$

END