

# 物理实验数学中心

Physics Expeiment Center




#### Download link: https://github.com/bliseu/phylab/

- [1]Qualitative study of atomic spectra, ROOM 532;
- [2] Equal thickness interference, ROOM 533;
- [3]Franck-Hertz experiment, ROOM 427;
- [4]Sound Velocity Measurement, ROOM 418;
- [5]11-cord potential difference meter, ROOM 417;
- [6]Amperemeter modification, ROOM 416;
- [7] Dielectric constant measurement, ROOM 429.

# Qualitative study of atomic spectra

Li Bin NJUPT

#### I. Diffraction grating



**Grating equation:** 

 $d*\sin\theta = k*\lambda, k=0, \pm 1, \pm 2....$ 

d: grating constant


θ: diffraction angle

k: order


**λ:** wavelength

# II. The adjustment of spectrometer


- >Final vision: see right
- >Method:
  - Three adjusting screw button under the loading platform
  - The telescope elevation adjusting screw
- >Steps:
  - Coarse adjustment
  - Fine ajustment





# collimator



#### Measurement of the first-order diffraction angles








Right side

Left side

## Read the angle



#### III. Experiment contents

1. With grating as diffraction device, measure the first-order diffraction angle of helium, determine the diffraction angle  $(\phi)$ -wavelength( $\lambda$ ) relation graph according to the known helium spectral wavelength:

Color	Wavelength (nm)	Color	Wavelength (nm)	Color	Wavelength (nm)		
Red(da rk)	706.52	Green(lig ht)	504.77	Blue	471.31		
Red	667.82	Green	501.57	Purple	447.15		
Yellow	587.56	Cyan	492.19	Purple (dark)	438.79		

at least 6 groups of data (uncovered rows), plot  $\lambda$ - $\varphi$  relation graph

Color	$\lambda(nm)$	$oldsymbol{ heta_{\!\scriptscriptstyle M}}$	$oldsymbol{ heta_{\!N}}$	$oldsymbol{ heta_{\!M}'}$	$oxed{ heta_{\!N}'}$	$\varphi = \frac{1}{4} \left( \left  \theta_{\scriptscriptstyle M}' - \theta_{\scriptscriptstyle M} \right  + \left  \theta_{\scriptscriptstyle N}' - \theta_{\scriptscriptstyle N} \right  \right)$
Red(d ark)	706.52					
Red	667.82					
Yellow	587.56					
Green( light)	504.77					
Green	501.57					
Cyan	492.19					
Blue	471.31					
Purple	447.15					
Purple (dark)	438.79					

2. Measure the diffraction angle of mercury lamp, find out the corresponding wavelength( $\lambda$ ) of mercury according to the former  $\lambda$ - $\varphi$  relation graph from the helium data.

Color	$ heta_{\!\scriptscriptstyle M}$	$ hilde{ heta_{\!\scriptscriptstyle \Lambda}}$	$\theta_{M}^{\prime}$	$\theta_N'$	$\varphi$ :	] = - 	$\frac{1}{1}\left( \left   heta  ight.  ight.  ight.$	<i>M</i>	$\left   heta_{\scriptscriptstyle M}  ight $	+	$\Theta_N'$ -	$-\left  heta_{\scriptscriptstyle N} ight $	)	$\lambda(nm)$
Yellow														
Green														
Blue													STATE OF	
Purple											N. S.			

## V. Original data

### 1. Helium lamp (six groups of data (uncovered)):

Color	$\lambda(nm)$	$\theta_{\!\scriptscriptstyle M}$	$oldsymbol{ heta_{\!N}}$	$\theta_{\!\scriptscriptstyle M}^{\prime}$	$  heta_N' $	$arphi = rac{1}{4} \left( \left   heta_{\scriptscriptstyle M}' -  heta_{\scriptscriptstyle M}  ight  + \left   heta_{\scriptscriptstyle N}' -  heta_{\scriptscriptstyle N}  ight   ight)$
Red(dark)						
Red	667.82			33		
Yellow	587.56			300		
Green(lig ht)						
Green	501.57					
Cyan	492.19					
Blue	471.31					
Purple	447.15					
Purple (dark)	438.79					

#### 2. Mercury lamp (four groups of data):

Color	$oxed{ heta_{\!M}}$	$oxed{ heta_{\!\scriptscriptstyle N}}$	$ heta_{\!\scriptscriptstyle M}'$	$oxed{ heta_{\!N}'}$	$\varphi$	) =	$\frac{1}{4}\left( \left  \cdot \right  \right)$	$ heta'_{\scriptscriptstyle M}$ -	- $\left.  heta_{\scriptscriptstyle M}  ight $	+	$ heta_{\scriptscriptstyle N}'$	$-\theta_{\scriptscriptstyle N}$	v )	$\lambda(nm)$
Yellow														
Green														
Blue														
Purple													À	

# END