

物理实验教学中心

Physics Experiment Center




Spectrometer

Li Bin
NJUPT

Experimental purposes

- 1、 The structure and adjustment of spectrometer.
- 2、 The measurement of the vertex angle of prism.

Instruments



Telescope

platform


collimator

Adjustments of spectrometer

Make the telescope focused at infinity




Make the axis of telescope & platform are perpendicular to the center axis



Make the axis of collimator are perpendicular to the center axis


Telescope

❖ Eyepiece (eye lens) focusing




❖ Objective lens focusing:

Make the double-faced mirror attached to the objective lens, move the eyepiece sleeve, then tighten the locking screw.




Vertical adjustment


- Final vision: see right
- Method:
 - ❖ Three adjusting screw button under the loading platform
 - ❖ The telescope elevation adjusting screw
- Steps:
 - ❖ Coarse adjustment
 - ❖ Fine adjustment




collimator



①




②



Measure the vertex angle

$$A = \frac{\varphi}{2} = \frac{1}{4} [|\theta'_M - \theta_M| + |\theta'_N - \theta_N|]$$




Calculate the vertex angle :

$$A = \frac{1}{4} (|\theta'_M - \theta_M| + |\theta'_N - \theta_N|)$$

We need to measure four angles:

$$\theta_M \quad \theta'_M \quad \theta_N \quad \theta'_N$$

Reading system



Data

➤Table

Angle NO.	θ_M	θ_N	θ_M'	θ_N'	$ \theta_M - \theta_M' $	$ \theta_N - \theta_N' $	$A = [\theta_M - \theta_M' + \theta_N - \theta_N']/4$	Average value \bar{A}
1	319°11'	139°8'	199°4'	19°4'				
2	54°39'	234°39'	174°53'	354°53'				
3	354°45'	174°44'	474°51'	294°47'				

END