

# 物理实验教学中心

*Physics Experiment Center*



# Dielectric constant measurement

Li Bin

NJUPT


## Purposes:

- 1、 Handle the principles and method of the measurement of dielectric constants.**
- 2、 Learn the method of data processing.**

## Principles:


### 1、 Measure the relative dielectric constant of solid dielectrics by alternating current bridge

In the air



$$C_1 = C_0 + C_{B1} + C_{F1}$$

Insert dielectric sample



$$C_2 = C_C + C_{B2} + C_{F2}$$

We have these functions:

$$C_o = \frac{\varepsilon_o S}{D}$$

$$C_c = \frac{\varepsilon_r \varepsilon_o S}{t + \varepsilon_r (D-t)}$$

$$C_{B1} = C_{B2}$$

$$C_{F1} = C_{F2}$$

$$C_c = C_2 - C_1 + C_o$$

$$\varepsilon_r = \frac{C_c \cdot t}{\varepsilon_o S - C_c (D-t)}$$

## 2、 Calculations of air dielectric constant and capacitance:

Vacuum dielectric constant  $\epsilon_0$ ,  $S_0$  is plate area,  $D$  is distance between plates, the capacitance  $C$  is shown as:

$$C = \frac{\epsilon_0 S_0}{D} + C_F$$

Set  $C = y$ ,  $\frac{1}{D} = x$


$$y = a + bx$$

We got  $a = C_F$ ,  $b = \epsilon_0 S_0$

Finally, we have  $\epsilon_0 = b/S_0$

## Instruments:

### 1、 Solid electrode



## 2、 Alternating current bridge

We choose 10kHz  
and read from  
'DISPLAY A' section.





### 3、 Solid dielectric sample


Solid samples are disk-shaped PTFE plastics, the material with high temperature resistance, corrosion resistance, high dielectric strength characteristics.



# Steps and Contents:

## 1、 Measure the dielectric constant of solid dielectrics by alternating current bridge

- (a). Measure the diameters of the sample as  $d$ , three times using Vernier caliper. Measure the thickness of sample as  $t$  three times using micrometer.
- (b). Let the two plates touch each other, and read the starting value of the micrometer as  $D_0$ .
- (c). Set the distance of the plates as 5.000mm. Note that the real value in micrometer should be  $5.000\text{mm} + D_0$ . Read the capacitance data from 'DISPLAY A' section as C1.
- (d). Insert the sample carefully, read the capacitance again as C2.
- (e). Repeat steps (d) and (e) twice.



## 2、 Measure the vacuum dielectric constant and capacitance

Remove the sample, set the plates distance to **1.000、  
1.100 ... 1.900 mm, and**

Read the corresponding capacitance  $C$ .

$$C = \frac{\epsilon_0 S_0}{D} + C_F$$

Set  $x = 1/D$ ,  $y = C$

Calculate  $C_F$  and  $\epsilon_0$  using origin software.

# Data processing:

## Table I

Starting value  $D_0 = \underline{0.000}$  mm, Distance  $D = \underline{5.000}$  mm

$\epsilon_0 = 8.85 \times 10^{-12} \text{ F/m}$

	D/mm	d/mm	t/mm	S/mm <sup>2</sup>	C <sub>0</sub> /pF	C <sub>1</sub> /pF	C <sub>2</sub> /pF	C <sub>C</sub> /pF	$\epsilon_r$
1	5.000			$\pi(d/2)^2 =$					
2									
3									
Average	-								

**d:** diameter of sample,

**t:** thickness of sample,

$S = \pi(d/2)^2$

$$C_0 = \frac{\epsilon_0 S}{D}$$

$$C_C = C_2 - C_1 + C_0$$

$$\epsilon_r = \frac{C_C \cdot t}{\epsilon_0 S - C_C (D - t)}$$

1 F =  $10^{12}$  pF,

1 m =  $10^3$  mm

## Table II

$D/\text{mm}$	1.000	1.100	1.200	1.300	1.400	1.500	1.600	1.700	1.800	1.900
$1/D(\text{mm}^{-1})$										
$C/\text{pF}$										


$$C = \frac{\epsilon_0 S_0}{D} + C_F$$

Set  $x = 1/D$ ,  $y = C$ ,  
calculate  $\epsilon_0$ .

The Slope  $K = \epsilon_0 S_0$ ,  
 $\epsilon_0 = K / S_0$

The Intercept is  $C_F$ .

Note:  $S_0 = 21.61 \text{ cm}^2$ .



Here is the weblink to download the slide:

<https://github.com/bliseu/phylab/>

1. Please finish the table I and II in the slide.
2. Plot a  $1/D$ - $C$  line, determine the slope ( $K = \epsilon_0 S_0$ ) and intercept ( $C_F$ ) of the line, then calculate  $\epsilon_0$  and compare it with the theoretical value ( $\epsilon_0 = 8.85 \times 10^{-12} \text{F/m}$ ).
3. Complete the report.

The DEADLINE is April 4, 2024.

END