

物理实验教学中心

Physics Experiment Center



Website: <https://sourceforge.net/projects/phy-njupt/files/2017Spring/>

Oscilloscope


Li Bin
NJUPT

Experimental Goals

1. Adjustment and use of oscilloscope.
2. Learn to use oscilloscope to observe voltage waveform.
3. Observation of Lissajous figures.

GDS-2062 Oscilloscope





ON/STBY

Acquire

Display

Utility

Program

Cursor

Measure

Save/Recall

Auto test/Stop

Auto Set

Help

Run/Stop


Hardcopy



VARIABLE

F Menu

**Chanel
switch**



HORIZONTAL

◀ POSITION ▶



HORI
MENU

TIME/DIV



TRIGGER

LEVEL




MENU



F05A Signal generator



Frequency

Amplitude



Reset :  + 8

Chanel A:  + 


Chanel B:  + 



Contents and Steps:


1. Settings of Signal generator:

(1) Set the frequency :




+ value + unit

(2) Set the voltage:




+ value + unit


(3) Set the phase difference between channel A & B :



+





value +



2. Observe voltage waveform on Oscilloscope

(1) Power on;

(2) Open the channels (corresponding lights on), press  or adjust the **VOLTS/DIV** & **TIME/DIV** buttons, observe the voltage waveform.



Sensitivity of Y
axis (voltage)


Sensitivity of
x axis (time)

3. Lissajous figures

(1) Turn on two channels;

(2)  + **F5** On oscilloscope

Switch to channel B:




+




On signal generator

Switch to channel A:



+




Phase differences:

On signal generator



$N_y=6$; Lissajous figures

$N_x:N_y=2:3$;
 $F_x:F_y=N_y:N_x=3:2$



$N_x=4$;

TABLEs

$$U_p = \frac{1}{2} U_{p-p} = \sqrt{2} U$$

1. The voltage of sinusoidal signals

NO. \	Voltage (V)	Sensitivity of Y axis (V)	D _y	U _{p-p} (v)	U _p (v)	U (v)
1	2.5					
2	3.0					
3	4.5					
4	5.0					

2. The period of sinusoidal signals

NO. \	$f \text{ (Hz)}$	Sensitivity of X axis (ms)	D_x	$T \text{ (ms)}$
1	400			
2	2000			
3	8000			
4	15000			

3. Lissajous figures

$$f_x = 3 \text{ kHz}, f_y = 6 \text{ kHz}$$

Phase diff.	0°	90°	180°	270°	360°
Lissajous figures					

END