

物理实验教学中心

Physics Experiment Center



Spectrometer


Li Bin

NJUPT

Experimental purposes


- 1、 Master the structure and adjustment of spectrometer.
- 2、 Master the determination of refractive index of prism with the method of minimum deviation angle.

Instruments



Adjustments of spectrometer

Make the **telescope** focused at infinity




Make the axis of **telescope** & **platform** are perpendicular to the **center axis**



Make the axis of **collimator** are perpendicular to the center axis


Telescope

❖ Eyepiece (eye lens) focusing:




❖ Objective lens focusing:

Make the double-faced mirror attached to the objective lens, move the eyepiece sleeve, then tighten the locking screw.




Vertical adjustment


- Final vision: see right
- Method:
 - ❖ Three adjusting screw button under the loading platform
 - ❖ The telescope elevation adjusting screw
- Steps:
 - ❖ Coarse adjustment
 - ❖ Fine adjustment




collimator




①



②



Minimum deviation angle




Calculate minimum deviation angle :

$$\delta_{\min} = \frac{1}{4} [|\theta'_M - \theta_M| + |\theta'_N - \theta_N|]$$

So we just need to measure four angles:

$$\theta_M \quad \theta'_M \quad \theta_N \quad \theta'_N$$

Reading system



Refractive index of prism

$$n = \frac{\sin \frac{\delta_{\min} + \alpha}{2}}{\sin \frac{\alpha}{2}}$$

$$\alpha = 60^\circ$$

END