

物理实验教学中心

Physics Experiment Center



Torsional pendulum


Li Bin
NJUPT

Experimental Purposes:

1. Determination of the **torsional constant** K of spiral spring using torsion pendulum;
2. Determination of the **rotational inertia** of different objects with a torsion pendulum, and compared with the theoretical value;
3. Verification of the **parallel axis theorem**.



Structure of torsional pendulum



Principles

1. Rotational inertia: I

2. Torsional pendulum

When the object is turned in a horizontal plane at a certain angle θ , the spring restoring moment M :

$M = -K*\theta$ K is the torsional constant .

And we also know that : $M = I * \beta$,

where $\beta = d^2\theta/dt^2 = -K*\theta/I = -\omega^2*\theta$,


the vibration period $T = 2\pi/\omega = 2\pi(I/K)^{1/2}$

3. Parallel axis theorem

$$I_c = I_0 + mx^2$$



Instruments:



A hollow metal cylinder, a solid plastic cylinder, a metal rod, and two small metal cylinders

Instrument panel



Adjustments :


- 1. Horizontal adjustment using air-bubble level;**
- 2. Light blocking rod can block the infrared emitting and receiving hole of the photoelectric detector.**
- 3. The swing angle should be between 40 - 90 degrees.**

Steps:


1. Determine the vibration periods ($10T$) of the loading plate and with cylindrical, calculate the torsion constant K ;



2. Determine the vibration periods ($10T$) of the hollow metal cylinder and the mental rod.



3. Verification of the parallel axis theorem;



4. Other measurements:
Electronic balance-----mass,
Tape-----length of the metal rod,
Vernier caliper-----Inside and outside diameters.

Post-processing:

1. Calculate K;

2. Calculate the experimental and theoretical rotational inertia I_{exp} and I_{theo} ,

3. Calculate the percentage error η :

$$\eta = |I_{\text{exp}} - I_{\text{theo}}| / I_{\text{theo}} * 100\%$$

Table I

$$k = 4\pi^2 \frac{I_1'}{\bar{T}_1^2 - \bar{T}_0^2} = \underline{3.699 \times 10^{-2} N/m},$$

$$\frac{K}{4\pi^2} = \underline{9.3685 \times 10^{-4} N/m}$$

Percentage error = $|I_{\text{exp}} - I_{\text{theo}}| / I_{\text{theo}} \times 100\%$

Items	m (kg)	Geometry (mm)		Period (s)		I _{theo} (kgm ²)	I _{exp} (kgm ²)	Percentage Error
Empty Plate	/	/		10T ₀	7.397	/	$I_0 = \frac{I_1'\bar{T}_0^2}{T_1^2 - T_0^2}$ 5.1288 × 10 ⁻⁴	/
					7.398			
					7.402			
				\bar{T}_0	0.739 9			
Plastic Cylinder	0.715	D ₁	100.06	10T ₁	12.26	$I_1' = \frac{1}{8}m\bar{D}_1^2$ =8.9459 × 10 ⁻⁴	$I_1 = \frac{K\bar{T}_1^2}{4\pi^2} - I_0$ =8.9449 × 10 ⁻⁴	0.012 %
			100.04		12.25			
			100.04		12.25			
		\bar{D}_1	100.04 7	\bar{T}_1	1.225 7			
Hollow Cylinder	0.701	D _{out}	100.10	10T ₂	15.17	$I_2' = \frac{1}{8}m(\bar{D}_{out}^2 + \bar{D}_{in}^2)$ =??	$I_2 = \frac{K\bar{T}_2^2}{4\pi^2} - I_0$ =??	??
			100.10		15.16			
			100.12					
		\bar{D}_{out}	100.10 7		15.18			
			D _{in}	93.73				
		93.74						
		93.75		\bar{T}_2				
		\bar{D}_{in}	93.74					
Metal Rod	0.134	L	609	10T ₄	21.06	$I_4' = \frac{1}{12}mL^2$ =??	$I_4 = \frac{K}{4\pi^2}\bar{T}_4^2$ =??	??
					21.06			
					21.05			
				\bar{T}_4	21.05 7			

$$k = 4\pi^2 \frac{I_1'}{\bar{T}_1^2 - \bar{T}_0^2} = \underline{3.699 \times 10^{-2} N/m},$$

$$\frac{K}{4\pi^2} = \underline{9.3685 \times 10^{-4} N/m}$$

$$m = \underline{238} \text{ g}$$

$$D_{\text{out}} = \underline{35.00} \text{ mm},$$

$$D_{\text{in}} = \underline{6.06} \text{ mm} \quad L_1 = \underline{33.04} \text{ mm}$$

$$I_5' = 2 \left[\frac{1}{16} m \left(D_{\text{out}}^2 + D_{\text{in}}^2 \right) + \frac{1}{12} m L_1^2 \right]$$

$$\text{Percentage error} = |I_{\text{exp}} - I_{\text{theo}}| / I_{\text{theo}} * 100\%$$

Table II

X (10 ⁻² m)	5. 00	10. 00	15. 00	20. 00	25. 00
10T(s)	24.06	31.03	40.13	49.97	60.22
	24.08	31.01	40.13	49.99	60.18
	24.07	30.08	40.10	49.95	60.20
\bar{T} (s)	2.407	3.101	4.012	4.997	6.020
I _{exp} (10 ⁻² kgm) $I = \frac{K}{4\pi^2} \bar{T}^2$	0.5422	??	??	??	??
I _{theo} (10 ⁻² kgm) $I' = I_4' + 2mx^2 + I_5'$	0.5413	0.8982	1.4932	2.3262	3.3972
Percentage Error	0.18%	??	??	??	??

END