

# 物理实验教学中心

*Physics Experiment Center*



# Michelson Interferometer

Li Bin  
NJUPT

- ◇ **Experimental Goals**
- ◇ **Experimental Principles**
- ◇ **Contents and Steps**
- ◇ **Data Processing**




# Experimental Goals


- ◇ Adjustment and application methods of the Michelson's Interferometer
- ◇ The principles and structure of the Michelson's Interferometer
- ◇ Observe the interference of equal inclination
- ◇ Measure the wavelength of red light from He-Ne laser

# Experimental Principles

## 1. The light path



The Michelson Interferometer  
light route



Interference of equal thick thin-film


## 2. The interference of equal inclination

The light path difference:

$$\delta = 2d \cos i' = \begin{cases} K\lambda & (K=0,1,2,\dots) \quad \text{Bright stripes} \\ (2K+1)\frac{\lambda}{2} & (K=0,1,2,\dots) \quad \text{Dark stripes} \end{cases}$$

(a) A greater thickness leads to denser stripes.

(b) When  $d$  increases, the stripes pop up; when  $d$  reduces, the stripes shrink in.



Regarding to central level of stripes (**i=0**) :

$$\left. \begin{array}{l} 2d = k\lambda \\ d = k\lambda/2 \\ d' = (k+1)\lambda/2 \end{array} \right\} \Rightarrow \Delta d = \Delta N \frac{\lambda}{2}$$

Wavelength:  $\lambda = \frac{2\Delta d}{\Delta N}$

# Contents and Steps

## 1. The condition to obtain the interference of equal inclination:


Two brightest light dots must be overlapped.

## 2. We observe the stripe diversification of equal inclination by changing $d$ .

( Stripes shrink or expand by rotating the coarse knob and precise knob.)




The  $M_1$  position is adjusted by coarse knob to remove M1 fleetly and precise knob to remove M1 slightly






3. To adjust circle center ( Adjust screws behind mirror M2 including horizontal and vertical screws )




4. Make the stripe clear and stable by adjusting coarse knob.

*Note: Before the measurement, M1's position was adjusted to be at about 50 mm.*



5. Count the circle number and read the data when the number reaches **50**.



Main ruler

**Attention:** The precise knob has to be rotated toward a fixed direction during measurement.

**Record:** There are five significant digits after the decimal point


Such as: **31.45 67 8mm**

Main ruler

Big knob

Small knob

Estimate



Reading window


Precise knob

# Operation video 1



迅捷视频剪辑软件

# Operation video 2



# Data Processing

## TABLE I

i	0	1	2	3	4	5
N	0	50	100	150	200	250
d <sub>i</sub> (mm)	50.07068	50.08669	50.10280	50.11896	50.13513	50.15144
Δd (mm)	Δd <sub>1</sub> =d <sub>3</sub> -d <sub>0</sub> =		Δd <sub>2</sub> =d <sub>4</sub> -d <sub>1</sub> =		Δd <sub>3</sub> =d <sub>5</sub> -d <sub>2</sub> =	
Average Δd (mm)	Δd=(Δd <sub>1</sub> +Δd <sub>2</sub> +Δd <sub>3</sub> )/3=					
λ=2Δd/ΔN (ΔN=150)						

Noted: wavelength( $\lambda$ ) of He-Ne laser is 632.8nm

END